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SGSOFT®

Wound dressing Instructions for use

Do not use if package is damaged

LOT

Batch number

Single use

STERILE R

Store at ambient temperature

Manufacturer

Used by Month and Year

0537

4

Manufacturer: Silvergreen Oy Teknobulevardi 3-5, 01530 Vantaa, Finland Issued: 21.8.2014

Intended use

The SGSOFT® wound dressing is intended for providing external mechanical protection to wounds and for absorbing exudates. The wound dressing also offers antimicrobial properties.

Product description

SGSOFT® is a sterile, disposable, non-woven wound dressing, which consists of viscose (VI) and polyester (PET) fibres containing ionic silver (Ag+). The ionic silver in the SGSOFT® wound dressing prevents bacterial colonisation in the wound dressing, thus preventing growth of bacteria in the wound.

Indications

The SGSOFT® wound dressing is for external use only. The wound dressing suits all exuding, superficial and deep wounds, such as slow to heal wounds, infected wounds or varicose ulcers.

Contraindications

The SGSOFT® wound dressing should not be used to treat internal wounds or on individuals who are known to have sensitivity to silver.

Instructions for use

1. Clean the wound in accordance with standard procedures.
2. Select a wound dressing which is slightly larger than the wound.
3. Open the sterile pack of wound dressings and fold or cut a wound dressing to a suitable size.
 - Place the wound dressing directly on the wound bed.
 - Secure the wound dressing and cover it with a secondary bandage, if needed.
 - Monitor the volume of exudates in the wound dressing and change the dressing, when needed.
4. If the wound dressing adheres to a wound producing little exudate, the dressing should be moistened with sterile saline solution prior to removing it from the wound.
5. Changing or removing the wound dressing.
 - Change the wound dressing when the volume of exudates exceeds the absorbency of the wound dressing.
 - Change or remove the wound dressing in accordance with the instructions given by a physician or healthcare personnel.
6. Dispose of the wound dressing in accordance with local waste management guidelines.

If the wound exudate volume is minimal, the wound dressing may be moistened with sterile saline solution or water. The interval of changing the wound dressing depends on the volume of exudate or the care instructions provided by a physician or healthcare personnel, but the dressing should be changed no later than after 7 days of use. The SGSOFT® wound dressing can be used under a compress.

Warnings

Care should be taken to not contaminate the wound dressing when opening the package. If the packaging material is damaged or opened or the wound dressing has been contaminated after opening the package, the wound dressing should not be used. The wound dressing must not be resterilised or reused in order to avoid a risk of infection.

Please note the product expiration date marked on the product packaging. Do not use an expired product.

When using wound dressings, wound healing should be monitored, in particular upon changing the wound dressing. In addition, the individual's overall condition and health should be monitored. If needed, a physician should be contacted.

Storage

The product should be stored and kept in a dry place.

Miscellaneous

The product does not contain any animal-derived substances or materials. The product does not contain latex.

SGSOFT® is a registered trademark
Instructions for Use IFU010108001 EN

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